

10X EIA Coating Buffer

Immunoassay plate coating solution

EIA coating buffer can be used with any EIA/ELISA system, including chemiluminescent and colorimetric. Provided as a 10X concentrate.

Advantages

- **RELIABLE** – quality-controlled coating buffer for reproducible results
- **VERSATILE** – compatible with chemiluminescent and colorimetric EIA/ELISA assays
- **CONVENIENT** – save time using pre-made buffer. No more weighing, mixing or pH adjusting.



Ordering Information

Catalog Number	Product	Size
R-03730-D25	10X EIA Coating Buffer	250 ml

Related Product

R-03728-E10	AdvanBlock™-EIA Blocking Solution	1 L
R-03024-D50	AdvanWash™ 10X Washing Solution	500 ml
K-16025-D25	ELISABright™, sufficient for twenty-five (25) 96-well ELISA plates	250 ml

Advansta Corporation

1505 Adams Drive, Suite B1 | Menlo Park, CA 94025

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com

Product information: www.advansta.com/products/EIA_Coating_Buffer

Copyright © 2016 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. AdvanBlock™, AdvanWash™ and ELISABright™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

www.advansta.com

